




# Implant

S Y S T E M

PRODUCT CATALOGUE



# INDEX.

<b>Advan implant solutions</b>	<b>4</b>
Excellence and innovation in implant dentistry	
<hr/>	
<b>ONE implant line</b>	<b>7</b>
ONE Internal	14
ONE Conical	26
<hr/>	
<b>Surgical Tools and Kits</b>	<b>39</b>
Precision and high performance in daily practice	
<hr/>	
<b>Implant line packaging</b>	<b>51</b>
Instructions for use, labels and implant cards	
<hr/>	
<b>Magellan guided Surgery</b>	<b>61</b>
From virtual planning to surgery	
<hr/>	
<b>Bibliography</b>	<b>68</b>

## ADVAN IMPLANT SOLUTIONS

---

The real value of an **implant system** becomes evident over time.

For over **30 years**, ADVAN implant solutions have been synonymous with quality, effectiveness and long-term aesthetic stability.

Numerous publications and extensive long-term clinical experience have shown that ADVAN maintains the stability of hard and soft tissues and ensures excellent aesthetic results even many years after implant placement. <sup>[6][11][12]</sup>


**An all-Italian story** written with passion, dedication and a relentless pursuit of perfection.

**Ing. Mario Zearo**  
CEO Advan Srl



## **From Titanium Bar to Therapeutic Solutions**

Advan products arise from the contribution of medical and engineering know-how, creating the basis for innovative developments. Thanks to its ability to innovate and respond promptly to the growing needs of dentists, Advan has developed a product line characterised by outstanding quality, durability and the competitiveness that only absolute excellence can guarantee.



# Implant line ONE.

---

# IMPLANT SYSTEM ONE.

ONE stands for simplicity, compatibility and versatility.  
**The implant line accessible to all.**

Advan ONE is the easy, effective and state-of-the-art implant system where the study of details translates into maximum versatility of use. Thanks to its unique dual loop, the fixture is easy to position. The design of the ONE implant has been studied and designed to achieve high primary stability with both compact and soft bone, even in the case of immediate post-extraction implants.



## A unique implant design

Modern, easy-to-use implant screws engineered for excellent primary stability together with tissue maintenance.



## Native digital system

ONE is seamlessly integrated into the most widely used digital implant planning and CAD/CAM work environments.



The One line meets the needs of the most experienced professionals and also of those new to implant surgery, offering high **quality**, maximum **compatibility** and **advanced design**.



ONE, is the first **modular implant system** designed to simplify the work of dental practices.

## YO**UNIVERSAL**IMPLANTOLOGYSYSTEM

**Advanced solution suitable for all bone types.**


- 1 **Ease of** surgical approach
- 2 **Two connections:** conical and internal
- 3 **A single** external morphology
- 4 **A single** surgical kit
- 5 Process **optimisation**
- 6 **Reduction** of costs

# IMPLANT DESIGN: One Internal and One Conical.

Two connections,  
a single *implant design*.

The ONE system features a single surgical kit and a single external endosseous screw morphology, with a choice of two different connections: One Internal and One Conical.

Regardless of the connection type chosen, the prosthetic connection guarantees native platform switching on all implant diameters. In addition, the single connection for all implant diameters makes for streamlined and foolproof implant prosthetic steps.



## One Internal

Internal hexagonal connection, to easily handle most implant cases from single edentulism to total rehabilitation (connection compatible with the world's most popular internal connection system).

CONNECTION



## One Conical

Conical connection, to better manage peri-implant hard and soft tissues in aesthetic areas or when tissue availability requires accurate management of the gingival connective tissue (connection compatible with the Advan GTB system).

CONNECTION

## The implant line suitable for all workflows: *traditional and digital.*

Advan's long-standing expertise in the management of state-of-the-art digital workflows enables all professionals to upgrade their traditional clinic and laboratory workflow, or to request a seamless integration of the implant system into their existing workflow.

The ONE implant line combines perfectly with guided surgery, allowing the implantologist to operate with a minimal degree of invasion, with the certainty of the best implant placement


relative to the amount of available bone and respect for the delicate surrounding anatomical structures. Today, virtual dental implant planning enables a prosthetically guided approach that can lead to the best possible prosthetic result in terms of design, aesthetics, occlusion and load.

**ONE prosthetics is suitable for even the most innovative digital workflows!**



### Scan

Simplify your workflow and plan the surgical interventions.



### Plan

Libraries available for all major CAD/CAM planning and modelling software.



### Treat

Safe treatments on customised projects.



### Monitor

Treatment monitoring and patient care.

## The digital workflow lets you:


- 1 Increase patient numbers by optimising dental procedures
- 2 Reduce stress with more precise planning, thanks to predictable workflows
- 3 Get the expected results by following the most direct and simplified workflow
- 4 Create precise, easy-to-clean restorations with long-lasting aesthetic results

MAIN FEATURES

### Cylindrical profile and machined platform

Design to optimise the load transmitted to the marginal bone.

- » Minimal cortical bone load <sup>[1][7][8]</sup>
- » Juxta-crestal positioning



### Fluid discharge channel

It distributes the blood clot over the entire implant surface, improving implant osseointegration.


- » Decreased hydrostatic compression of blood clotting and post-operative swelling
- » The clot, the best osteoinductive factor, in contact with the entire surface



### Surface treatment OsseoGRIP®

Non-aggressive treatment on the implant surface: sand-blasting with hydroxyapatite and removal by acidification with the exception of the prosthetic platform of the first coronal section for safe subperiosteal positioning.<sup>[3]</sup>

- » Increased BIC compared to other surface treatments, for faster bone regrowth <sup>[1][12]</sup>
- » Low surface roughness




### Double conical profile of the plant core

The cylindrical profile of the implant neck allows minimal load on the cortical bone.<sup>[1]</sup>

- » Ideal chewing load distribution
- » Suitable for immediate loading
- » High primary stability with possibility of soft bone sub-preparations





### Variable profile progressive threading

Double-lead threading with variable profile, aggressive in the apical section and trapezoidal in the coronal section.

- >> Suitable for bicortical protocols
- >> Suitable for immediate post-extraction



### Horizontal platform switching

The prosthetic connection, regardless of the type of connection chosen, provides native platform switching on all implant diameters, with the advantage of increased space available for peri-implant soft tissue.<sup>[6][11][12][13]</sup>

- >> Prosthetic mis-matching at the implant connection level
- >> Increased volume available for peri-implant soft tissues <sup>[6][11][12][13]</sup>



### Active design of the implant apex

Atraumatic convex apex and exceptional cutting capacity for greater stability during insertion.

- >> Reduced risk of damaging the more delicate anatomical structures
- >> Ease and speed of insertion

**The ONE dental implant** is a system with controlled and gradual transmission of the chewing load.

# One Internal CONNECTION.



One Internal is designed with the most common **internal hexagon connection** on the market. Simple and intuitive, it guarantees reliability and safety in the surgical act.



## Internal hexagon

Internal hexagon connection. Prevents rotation between implant and abutment. It provides an indication of positioning and ensures precise and firm insertion using dedicated tools.

## Prosthetic connection

A single prosthetic connection for all implant diameters, offering a wide range of prosthetic solutions.

## Platform Switching


Horizontal platform switching improves the distribution of the chewing load and increases the volume available for peri-implant soft tissue by stabilising the papilla.<sup>[6][11][12]</sup>

## OsseoGRIP®

Machined platform and OsseoGRIP® surface treatment down to 0.5 mm from the prosthetic platform for secure subperiosteal positioning.

# Size of the implant.

Neck diameter ▶



Apical core diameter ▶

Apical loop diameter ▶




	Ø 3.75	Ø 4.2	Ø 4.65
		<b>N4206</b> L 6 mm	<b>N4706</b> L 6 mm
	<b>N3808</b> L 8 mm	<b>N4208</b> L 8 mm	<b>N4708</b> L 8 mm
	<b>N3810</b> L 10 mm	<b>N4210</b> L 10 mm	<b>N4710</b> L 10 mm
	<b>N3811</b> L 11.5 mm	<b>N4211</b> L 11.5 mm	<b>N4711</b> L 11.5 mm
	<b>N3813</b> L 13 mm	<b>N4213</b> L 13 mm	<b>N4713</b> L 13 mm

The One internal implant screw is colour-coded blue for better recognisability than the conical connection in both the information materials and packaging (specification on page 51). For further information on the implant, packaging and surgical procedures please refer to the ONE internal Surgical Manual.


# Cover screw and healing abutments.

## Cover screw

The surgical cover screw is intended to be temporarily screwed into the implant connection in cases where primary stability is not achieved. Packaged sterile. Made of **Ti gr.23**.




Use the prosthetic driver.  
Tightening torque **10 Ncm**.




**01NVT02**  
Cover screw  
H 0.0 mm  
(3 pieces)

## Healing abutments

They allow proper conditioning of the mucosal tunnel. All profiles are of a larger emergent diameter compared to the corresponding prosthetic abutment to allow convenient and comfortable positioning for the patient without the need for anaesthetic. Packaged sterile. Made of **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **20 Ncm**.



**05NMG12**  
Healing abutment  
Ø 5.1 mm / H 2.0 mm

**05NMG13**  
Healing abutment  
Ø 5.1 mm / H 3.5 mm

**05NMG14**  
Healing abutment  
Ø 5.1 mm / H 4.5 mm



# Abutment Transfer, Analog and Digital Analog Abutments.

## Transfer Abutments


Impression abutments for use with open or closed technique:

- Open tray impression abutment with large retention portion with drains to facilitate ligation of several transfer abutments.
- Cannula for extending the fixation screw, for open technique: lets the user facilitate the unscrewing of the fixation screw in the event that the elastomer covers its head. Made of PMMA.
- Abutment for closed tray impression taking, specially shaped to maximise stability in the impression material and make both removal of the tray and subsequent repositioning of the abutment simple and safe.

The transfer abutment and implant attachment screw are made of **Ti gr.23**.




Use the prosthetic driver.  
Tightening torque **7 Ncm**.



**01INTO10**  
Open tray  
transfer abutment  
H 16.2 mm




**01CVT10**  
Double tray  
transfer cannula  
(6 pieces)



**01NTC10**  
Closed tray  
transfer abutment  
H 10.5 mm

## Traditional analog-digital

Analog replicating the implant connection. Component designed to also be used as a digital analog.  
Made of **Ti gr.23**.



**01NAN10**  
Analog

## InLab digital analog

The digital analog lets you work in the InLab software environment. Standard platform Ø 4,2 mm.  
Made of **Ti gr.23**.




**01SAN02**  
Digital analog


# Scan abutment and Ti-Base abutments.

## Scan abutment

Abutment for impression taking.  
Made of **sandblasted Ti gr.23**.




Use the prosthetic driver.  
Tightening torque **20 Ncm**.  
Prosthetic screw **TYPE 2 (01NFL18)** supplied.




## Ti-Base abutments

Ti-Base is the titanium base for CAD/CAM rehabilitation; it contributes to a cost-effective workflow for the production of customised abutments. Accurate digital acquisition of implant position by body scan: extraoral on the model or intraoral. Adhesive connection of Ti-Base with the mesostructure or abutment. Suitable for the Cerec workflow (Standard S) by Dentsply Sirona. Made of **Ti gr.23**.




Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 1 (01NFC18)** supplied.




# Uni-Base Abutments and Overcast Abutments with CoCr28Mo Base.

## Uni-Base abutments


Uni-Base is the titanium base for CAD/CAM rehabilitation; it contributes to a cost-effective workflow for the production of customised abutments. Accurate digital acquisition of implant position by body scan: extraoral on the model or intraoral. Adhesive connection of Uni-Base with the mesostructure or abutment. Made of **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 1 (01NFC18)** supplied.




**41NUB20**  
Uni-Base abutment  
H 0.8 mm / L 8.9 mm  
Locking




**41NUB30**  
Uni-Base abutment  
H 0.8 mm / L 8.9 mm  
Not Locking

## Castable abutments with CoCr28Mo base


Overcast abutment in castable acrylic with milled CoCr28Mo base. The abutment is designed for single crown or implant bridge rehabilitations and is commonly indicated for bar-supported overdentures. The precisely machined CoCr28Mo base offers an absolutely precise fit to the implant. The abutments can be individually milled and shaped for a better emergence profile and an aesthetically natural result. It is directly connected to the implant via the fastening screw. PMMA castable cylinder.



Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 1 (01NFC18)** supplied.



**41MCN40**  
CoCr28Mo base Castable abutment  
Ø 4.1 mm / H 0.8 mm / L 5.0 mm  
Locking




**41MCR40**  
CoCr28Mo base Castable abutment  
Ø 4.1 mm / H 0.8 mm / L 5.0 mm  
Not Locking




# Titanium abutments.

## Titanium abutments


Pre-finished abutments for rapid, highly aesthetic prostheses. The predefined angle of the abutment makes it easy to correct lack of implant parallelism. The abutment pin has a flat anti-rotational face and an anatomical profile to reduce bulk. Made of **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 1 (01NFC18)** supplied.




Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 1 (01NFC18)** supplied.  
Angled to the face of the hexagon.



# Multi Unit Abutment.

## MUA

Multi-prosthetic screw-retained connection system for screw-retained bridges, screw-retained restorative bars, Toronto-Bridge, rehabilitations of totally edentulous arches. Ideal for rehabilitations with more than 3 mm of mucous membrane. Sterile packaged with pre-assembled carrier. Made of **Ti gr.23**.



**05MND01**  
MUA - Straight  
Ø 4.8 mm / H 1.0 mm / α 0°  
Not Locking

**05MNA11**  
MUA - Angled  
Ø 4.8 mm / H 1.0 mm / α 17°  
Locking

**05MNA32**  
MUA - Angled  
Ø 4.8 mm / H 2.0 mm / α 30°  
Locking

**05MND02**  
MUA - Straight  
Ø 4.8 mm / H 2.0 mm / α 0°  
Not Locking


**05MNA12**  
MUA - Angled  
Ø 4.8 mm / H 2.0 mm / α 17°  
Locking

**05MNA33**  
MUA - Angled  
Ø 4.8 mm / H 3.5 mm / α 30°  
Locking


**05MND03**  
MUA - Straight  
Ø 4.8 mm / H 3.5 mm / α 0°  
Not Locking

**05MNA13**  
MUA - Angled  
Ø 4.8 mm / H 3.5 mm / α 17°  
Locking

**05MND04**  
MUA - Straight  
Ø 4.8 mm / H 4.5 mm / α 0°  
Not Locking



Use MUA drivers.  
Tightening torque **30 Ncm**.  
Pre-assembled carrier included.



Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 3 (01NFM18)** supplied.  
Angled to the face of the hexagon.  
Pre-assembled carrier included.




Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 4 (02NFM18)** supplied.  
Angled to the face of the hexagon.  
Pre-assembled carrier included.


# MUA prosthetic components.

## MUA Prosthetic Components


The MUA system features a series of prosthetic components that allow soft tissue conditioning, making an impression with an open technique, making a master model with a dedicated analogue and making provisional and definitive prosthetic rehabilitations.




Use the prosthetic driver.  
Tightening torque **7 Ncm**.  
Made of **Ti gr 23**.




**01TMU10**  
MUA - Transfer abutment  
H 11.0 mm




**01AMU10**  
MUA - Analog




**01SMU10**  
MUA - Scan abutment  
Prosthetic screw TYPE D  
(01-VE14) included  
H 8.0 mm




Use the prosthetic driver.  
Tightening torque **15 Ncm**.  
Prosthetic screw **TYPE D**  
(01-VE14) included.  
Made of **Ti gr 23**.




**01CMU20**  
MUA - Temporary  
Titanium sleeve  
H 12.0 mm




**01GMU10**  
MUA - Healing cap  
H 4.3 mm




**01BMU10**  
MUA - Uni-Base  
H 7.0 mm



**01BMU20**  
MUA - Uni-Base  
H 9.0 mm



Use the prosthetic driver.  
Tightening torque **15 Ncm**.  
Prosthetic screw **TYPE D**  
(01-VE14) included.  
Made of **PMMA**.




**01CMU10**  
MUA - Castable sleeve  
H 12.0 mm


# Pre-Milled Abutments and Prosthetic Screws.

## Pre-Milled Abutments

This prosthetic component has a pre-turned connection with all the features and strict tolerances of the ADVAN implant system. The Pre-Milled abutment makes it possible to obtain a secondary component with a customised transgingival design, regardless of the implant position. Made of **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **30 Ncm**.  
Prosthetic screw **TYPE 1 (01NFC18)** supplied.



## Prosthetic Screws

Screws for prosthetic components.  
Made of **Ti gr.23**.



**01NFC18**  
Prosthetic screw  
TYPE 1 M1.8  
(3 pieces)



**01NFL18**  
Prosthetic screw  
TYPE 2 M1.8  
(3 pieces)



**01NFM18**  
Prosthetic screw  
TYPE 3 M1.8  
(3 pieces)



**02NFM18**  
Prosthetic screw  
TYPE 4 M1.8  
(3 pieces)



**01-VE14**  
Prosthetic screw  
TYPE D M1.4  
(3 pieces)

## Dynamic Screw TYPE 1

Dynamic screw for prosthetic components to angle the prosthetic hole. Use only with dynamic driver. Dynamic drivers are identified by a yellow colour-coded label. Made of **Ti gr.23**.




**01-EE20**  
Dynamic screw  
TYPE 1 M1.8  
(3 pieces)

# Retentive Multi-prosthetic System OT EQUATOR®.

## Retentive Prosthetic Multisystem

**OT EQUATOR**

*Smart*  
**BOX**



### 50° DIVERGENCE CORRECTION

OT Equator retentive prosthetic multisystem, spare parts, accessories and specific tools are made by **RHEIN83**® srl



## Retentive Multi-prosthetic System OT Equator®

OT Equator® overdenture abutments provide excellent retention with minimal vertical and horizontal space requirements. The retention of OT Equator® attachments can be adapted to the needs of patients by choosing between different OT Equator® or NORMO retention caps. SmartBOX cap boxes allow correction of up to 50° of misalignment between implants, ensuring safe and comfortable prosthesis placement. Made of **Ti gr. 5** and coated with **TiN**.


### OT EQUATOR® SMARTBOX PACK

**Content:**

- 1 OT Equator® abutment
- 1 SmartBox container
- 1 Protective disc for fitting
- 1 Protective cap
- 4 Assorted retention caps  
(1 extra-soft, 1 soft, 1 standard, 1 strong)



Use OVD drivers.  
Tightening torque **35 Ncm**.



**131CV1**  
Ø 3.5 mm / H 1.0 mm

**131CV2**  
Ø 3.5 mm / H 2.0 mm

**131CV3**  
Ø 3.5 mm / H 3.0 mm

**131CV4**  
Ø 3.5 mm / H 4.0 mm

**131CV5**  
Ø 3.5 mm / H 5.0 mm



**330SBE**  
Cap holder



**140CEV**  
Retentive cap  
Purple Strong  
Kg 2.7  
(4 pieces)



**140CET**  
Retentive cap  
White Standard  
Kg 1.8  
(4 pieces)



**140CER**  
Retentive cap  
Pink Soft  
Kg 1.2  
(4 pieces)



**140CEG**  
Retentive cap  
Yellow Extra Soft  
Kg 0.6  
(4 pieces)



**140CEN**  
Retentive cap  
Black  
(for laboratory)  
(4 pieces)

# CONNECTION

## One Conical.



The One Conical implant line is designed with the **conical connection** of the GTB implant line, which guarantees the absence of micro movements.<sup>[2][4]</sup>



### Conical connection

One conical has the same sealing conical connection as Advan's GTB implant line. There are no detectable micro-gaps at the implant-abutment interface in the connection area (<1 µm).<sup>[2][4]</sup>

### OsseoGRIP®

Machined platform and OsseoGRIP® surface treatment down to 0.3 mm from the prosthetic platform for secure subperiosteal positioning.

### Single connection diameter

Ease of use thanks to the unique fixture design and the adoption of a single prosthetic platform for all implant diameters.

### Platform Switching


Horizontal platform switching improves the distribution of the chewing load and increases the volume available for peri-implant soft tissue by stabilising the papilla.<sup>[6][11][12][13]</sup>

### Octagonal anti-rotational index

The anti-rotational octagon prevents rotation between implant and abutment.

# Size of the implant.

Neck diameter ►



	Ø 3.3	Ø 3.75	Ø 4.2	Ø 4.65
			<b>C4206</b> L 6 mm	<b>C4706</b> L 6 mm
	<b>C3308</b> L 8 mm	<b>C3808</b> L 8 mm	<b>C4208</b> L 8 mm	<b>C4708</b> L 8 mm
	<b>C3310</b> L 10 mm	<b>C3810</b> L 10 mm	<b>C4210</b> L 10 mm	<b>C4710</b> L 10 mm
	<b>C3311</b> L 11.5 mm	<b>C3811</b> L 11.5 mm	<b>C4211</b> L 11.5 mm	<b>C4711</b> L 11.5 mm
	<b>C3313</b> L 13 mm	<b>C3813</b> L 13 mm	<b>C4213</b> L 13 mm	<b>C4713</b> L 13 mm

The One conical implant screw is colour-coded green for better recognisability compared to the internal connection both in the information materials and in the packaging (specification on page 51). For further information on the implant, packaging and surgical procedures refer to the ONE conical Surgical Manual.

# Cover screw and healing abutments.

## Cover screw

The surgical cover screw is intended to be temporarily screwed into the implant connection in cases where primary stability is not achieved. Packaged sterile. Made of **Ti gr.23**.




Use the prosthetic driver.  
Tightening torque **7 Ncm**.



**01GVT02**  
Cover screw  
H 0.0 mm  
(3 pieces)

## Healing abutments

They allow proper conditioning of the mucous canal. All profiles are of a larger emergent diameter compared to the corresponding prosthetic abutment to allow convenient and comfortable positioning for the patient without the need for anaesthetic. Packaged sterile. Made of **Ti gr.23**.




**04GMG02**  
Healing abutment  
Ø 4.0 mm / H 2.0 mm

**04GMG03**  
Healing abutment  
Ø 4.0 mm / H 3.5 mm

**04GMG04**  
Healing abutment  
Ø 4.0 mm / H 4.5 mm



Use the prosthetic driver.  
Tightening torque **7 Ncm**.



**050MG12**  
Healing abutment  
Ø 5.1 mm / H 2.0 mm

**050MG13**  
Healing abutment  
Ø 5.1 mm / H 3.5 mm

**050MG14**  
Healing abutment  
Ø 5.1 mm / H 4.5 mm


# Abutment Transfer, Analog and Digital Analog.

## Transfer Abutments

Impression abutments for use with open or closed technique:

- Open tray impression abutment with large retention portion with drains to facilitate ligation of several transfer abutments.
- Optional transfer screws available in two different heights.
- Abutment for closed tray impression taking, specially shaped to maximise stability in the impression material and make both removal of the tray and subsequent repositioning of the abutment simple and safe.

The transfer abutment and implant attachment screw are made of **Ti gr.23**.



## Traditional analog-digital

Analog replicating the implant connection. A component designed to also be used as a digital analog. Made of **Ti gr.23**.



**01GAN10**  
Analog

## InLab digital analog

The digital analog lets you work in the InLab software environment. Standard platform Ø 4.2 mm. Made of **Ti gr.23**.




**01SAN01**  
Digital analog


# Scan abutment and Ti-Base abutments.

## Scan abutment

Abutment for impression taking.  
Made of **sandblasted Ti gr.23**.




Use the prosthetic driver.  
Tightening torque **7 Ncm**.  
Prosthetic screw **TYPE A (01VFM16)** supplied.




**04GSB10**  
Scan abutment  
L 10.5 mm

## Ti-Base abutments

Ti-Base is the titanium base for CAD/CAM rehabilitation; it contributes to a cost-effective workflow for the production of customised abutments. Accurate digital acquisition of implant position by body scan: extraoral on the model or intraoral. Adhesive connection of Ti-Base with the mesostructure or abutment. Suitable for the Cerec workflow (Standard S) by Dentsply Sirona. Made of **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **25 Ncm**.  
Prosthetic screw **TYPE A (01VFM16)** supplied.



**41GLC00**  
Ti-Base abutment  
Ø 4.1 mm / H 0.7 mm / L 5.6 mm


**41GLC02**  
Ti-Base abutment  
Ø 4.1 mm / H 2.0 mm / L 6.7 mm

**41NLC30**  
Ti-Base abutment  
Ø 4.1 mm / H 3.5 mm / L 8.3 mm

# Uni-Base Abutments and Overcast Abutments with CoCr28Mo Base.

## Uni-Base abutments


Uni-Base is the titanium base for CAD/CAM rehabilitation; it contributes to a cost-effective workflow for the production of customised abutments. Accurate digital acquisition of implant position by body scan: extraoral on the model or intraoral. Adhesive connection of Uni-Base with the mesostructure or abutment. Made of **Ti gr.23**.



<p><b>41CUB20</b> Uni-Base abutments Ø 4.1 mm / H 0.7 mm / L 9.0 mm Locking</p>	<p><b>41CUB30</b> Uni-Base abutments Ø 4.1 mm / H 0.7 mm / L 9.0 mm Not Locking</p>	<p><b>41UNN02</b> Uni-Base abutments Ø 4.1 mm / H 2.0 mm / L 6.7 mm Locking</p>	<p><b>41UNR02</b> Uni-Base abutments Ø 4.1 mm / H 2.0 mm / L 6.7 mm Not Locking</p>
<p><b>25 Ncm</b></p> <p>Use the prosthetic driver. Tightening torque <b>25 Ncm</b>. Prosthetic screw <b>TYPE A (01VFM16)</b> supplied.</p>	<p><b>25 Ncm</b></p> <p>Use the prosthetic driver. Tightening torque <b>25 Ncm</b>. Prosthetic screw <b>TYPE B (01VFL16)</b> supplied.</p>	<p><b>42CUB20</b> Uni-Base abutments Ø 4.1 mm / H 3.5 mm / L 8.3 mm Locking</p>	<p><b>42CUB30</b> Uni-Base abutments Ø 4.1 mm / H 3.5 mm / L 8.3 mm Not Locking</p>
		<p><b>41UNN04</b> Uni-Base abutments Ø 4.1 mm / H 4.5 mm / L 9.3 mm Locking</p>	<p><b>41UNR04</b> Uni-Base abutments Ø 4.1 mm / H 4.5 mm / L 9.3 mm Not Locking</p>

## Castable abutments with CoCr28Mo base

Overcast abutment in castable acrylic with milled CoCr28Mo base. The abutment is designed for single crown or implant bridge rehabilitations and is commonly indicated for bar-supported overdentures. The precisely machined CoCr28Mo base offers an absolutely precise fit to the implant. The abutments can be individually milled and shaped for a better emergence profile and an aesthetically natural result. It is directly connected to the implant via the fastening screw. PMMA castable cylinder.



Use the prosthetic driver.  
Tightening torque **25 Ncm**.  
Prosthetic screw **TYPE A (01VFM16)** supplied.


**41CCN40**  
CoCr28Mo base Castable abutment  
Ø 4.1 mm / H 1.0 mm / L 5.0 mm  
Locking

**41CCR40**  
CoCr28Mo base Castable abutment  
Ø 4.1 mm / H 1.0 mm / L 5.0 mm  
Not Locking




# Titanium abutments.

## Titanium abutments

Predefined abutments for fast, highly aesthetic prostheses. The predefined angle of the abutment makes it easy to correct lack of implant parallelism. The abutment pin has a flat anti-rotational face and an anatomical profile to reduce bulk. Made in **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **25 Ncm**.  
Prosthetic screw **TYPE A (01VFM16)** supplied.  
Angled compared to the vertex of the octagon.






# Multi Unit Abutment.

## MUA

Multi-prosthetic screw-retained connection system for screw-retained bridges, screw-retained restorative bars, Toronto-Bridge, rehabilitations of totally edentulous arches. Ideal for rehabilitations with more than 3 mm of mucous membrane. Sterile packaged with pre-assembled carrier. Made of **Ti gr.23**.



### 05GND00

MUA - Straight  
 Ø 4.8 mm / H 1.0 mm / α 0°  
 Not Locking

### 05GNA10

MUA - Angled  
 Ø 4.8 mm / H 1.0 mm / α 17°  
 Locking

### 05GNA32

MUA - Angled  
 Ø 4.8 mm / H 2.0 mm / α 30°  
 Locking

### 05GND02

MUA - Straight  
 Ø 4.8 mm / H 2.0 mm / α 0°  
 Not Locking

### 05GNA12

MUA - Angled  
 Ø 4.8 mm / H 2.0 mm / α 17°  
 Locking

### 05GNA33

MUA - Angled  
 Ø 4.8 mm / H 3.5 mm / α 30°  
 Locking

### 05GND03

MUA - Straight  
 Ø 4.8 mm / H 3.5 mm / α 0°  
 Not Locking

### 05GNA13

MUA - Angled  
 Ø 4.8 mm / H 3.5 mm / α 17°  
 Locking

### 05GND04

MUA - Straight  
 Ø 4.8 mm / H 4.5 mm / α 0°  
 Not Locking



Use MUA drivers.  
 Tightening torque **25 Ncm**.  
 Pre-assembled carrier included.



Use the prosthetic driver.  
 Tightening torque **25 Ncm**.  
 Prosthetic screw **TYPE A (01VFM16)** supplied.  
 Angled to the face of the octagon.  
 Pre-assembled carrier included.




Use the prosthetic driver.  
 Tightening torque **25 Ncm**.  
 Prosthetic screw **TYPE C (01VGF16)** supplied.  
 Angled to the face of the octagon.  
 Pre-assembled carrier included.


# MUA prosthetic components.

## MUA Prosthetic Components


The MUA system features a series of prosthetic components that allow soft tissue conditioning, making an impression with an open technique, making a master model with a dedicated analogue and making provisional and definitive prosthetic rehabilitations.




Use the prosthetic driver.  
Tightening torque **7 Ncm**.  
Made of **Ti gr 23**.




**01TMU10**  
MUA - Transfer abutment  
H 11.0 mm




**01AMU10**  
MUA - Analog




**01SMU10**  
MUA - Scan abutment  
Prosthetic screw TYPE D  
(01-VE14) included  
H 8.0 mm




Use the prosthetic driver.  
Tightening torque **15 Ncm**.  
Prosthetic screw **TYPE D**  
(01-VE14) included.  
Made of **Ti gr 23**.




**01CMU20**  
MUA - Temporary  
Titanium sleeve  
H 12.0 mm




**01GMU10**  
MUA - Healing cap  
H 4.3 mm




**01BMU10**  
MUA - Uni-Base  
H 7.0 mm



**01BMU20**  
MUA - Uni-Base  
H 9.0 mm



Use the prosthetic driver.  
Tightening torque **15 Ncm**.  
Prosthetic screw **TYPE D**  
(01-VE14) included.  
Made of **PMMA**.




**01CMU10**  
MUA - Castable sleeve  
H 12.0 mm


# Pre-Milled Abutments and Prosthetic Screws.

## Pre-Milled Abutments

This prosthetic component has a pre-turned connection with all the features and strict tolerances of the ADVAN implant system. The Pre-Milled abutment makes it possible to obtain a secondary component with a customised transgingival design, regardless of the implant position. Made of **Ti gr.23**.



Use the prosthetic driver.  
Tightening torque **25 Ncm**.  
Prosthetic screw **TYPE A (01VFM16)** supplied.



## Prosthetic Screws

Screws for prosthetic components.  
Made of **Ti gr.23**.



**01VFM16**  
Prosthetic screw  
TYPE A M1.6  
(3 pieces)



**01VFL16**  
Prosthetic screw  
TYPE B M1.6  
(3 pieces)



**01VGF16**  
Prosthetic screw  
TYPE C M1.6  
(3 pieces)



**01-VE14**  
Prosthetic screw  
TYPE D M1.4  
(3 pieces)

## Dynamic Screws

Dynamic screws for prosthetic components to angle the prosthetic hole. Use only with dynamic driver. Dynamic drivers are identified by a yellow colour-coded label. Made of **Ti gr.23**.



**01EVC20**  
Dynamic screw  
TYPE A M1.6  
(3 pieces)




**01EVL20**  
Dynamic screw  
TYPE B M1.6  
(3 pieces)

# Retentive Multi-prosthetic System OT EQUATOR®.

## Retentive Prosthetic Multisystem

**OT EQUATOR**

*Smart*  
**BOX**



### 50° DIVERGENCE CORRECTION

OT Equator retentive prosthetic multisystem, spare parts, accessories and specific tools are made by **RHEIN83**® srl

## Retentive Multi-prosthetic System OT Equator®

OT Equator® overdenture abutments provide excellent retention with minimal vertical and horizontal space requirements. The retention of OT Equator® attachments can be adapted to the needs of patients by choosing between different OT Equator® or NORMO retention caps. SmartBOX cap boxes allow correction of up to 50° of misalignment between implants, ensuring safe and comfortable prosthesis placement. Made of **Ti gr. 5** and coated with **TiN**.


### OT EQUATOR® SMARTBOX PACK

**Content:**

- 1 OT Equator® abutment
- 1 SmartBox container
- 1 Protective disc for fitting
- 1 Protective cap
- 4 Assorted retention caps  
(1 extra-soft, 1 soft, 1 standard, 1 strong)



Use OVD drivers.  
Tightening torque **35 Ncm**.



**131GER41**

Ø 4.0 mm / H 1.0 mm

**131GER42**

Ø 4.0 mm / H 2.0 mm

**131GER43**

Ø 4.0 mm / H 3.5 mm

**131GER44**

Ø 4.0 mm / H 4.5 mm

**131GER45**

Ø 4.0 mm / H 5.5 mm



**330SBE**  
Cap holder



**140CEV**  
Retentive cap  
Purple Strong  
Kg 2.7  
(4 pieces)



**140CET**  
Retentive cap  
White Standard  
Kg 1.8  
(4 pieces)



**140CER**  
Retentive cap  
Pink Soft  
Kg 1.2  
(4 pieces)



**140CEG**  
Retentive cap  
Yellow Extra Soft  
Kg 0.6  
(4 pieces)



**140CEN**  
Retentive cap  
Black  
(for laboratory)  
(4 pieces)



**Surgical tools**  
and kits.

---

# SURGICAL TOOLS.

## We make innovative surgical tools to improve the daily practice of professionals

Reliable and easy to use, ADVAN surgical tools are your ideal partner for restorative dental procedures. We offer high-performance precision surgical tools designed to meet your patients' expectations and restore even the most difficult dental cases safely and without interruption.

The drivers and drill stop sets are made of stainless steel for durability, ensuring simplified, accurate and successful dental implant procedures.





# SURGICAL KITS

## Advan.

### Surgical Kits

ADVAN surgical kits are designed with all the tools required to carry out, in a logical and sequential manner, the steps that characterise an optimal surgical procedure. The arrangement and colour coding of the tools allows them to be quickly identified.

The kit consists of a single box containing all the surgical and prosthetic tools necessary for performing the osteotomy, placing the implant and its prosthesis. Including the surgical box made of Radel (non-toxic, autoclavable, metal-free), the surgical kit must be combined with the Implant Driver Set dedicated to the type of implant connection chosen: one conical or one internal. A kit organised to give the clinician maximum freedom. A simple, logged workflow for each ADVAN implant line.

The ultimate in functionality and practicality in a subdivided and perfectly organised box: all tools of the ONE line are housed in a single container, organised according to a logical and intuitive path.



**ONE KIT CODE: OKC10**

#### THE SURGICAL KIT INCLUDES:

CODE	PRODUCT
07-FI20	Pointed starter drill
07-FI10	Lindemann drill
07-FP10	Marking drill
01-MP12	Double diameter depth probe (2 pieces)
07FAB20	Drill $\varnothing$ 2.0 mm
07FAB28	Step drill $\varnothing$ 2.4 / 2.8 mm
07FAB32	Step drill $\varnothing$ 2.8 / 3.2 mm
07FAB36	Step drill $\varnothing$ 3.2 / 3.65 mm


CODE	PRODUCT
07FAB41	Step drill $\varnothing$ 3.65 / 4.1 mm
07FAB45	Step drill $\varnothing$ 4.1 / 4.5 mm
02EGM10	Manual prosthetic driver medium L 10.0 mm
02EGM15	Manual prosthetic driver long L 15.0 mm
07-EG10	Prosthetic handpiece driver medium L 10.0 mm
02-CT20	Torque ratchet 10-70 Ncm
02-AC50	Ratchet adapter
010LT10	ONE implant-sizing overlays

# DRILLS

## ONE implant line.


### Drill selection

The drills in the ONE implant line feature a single-ring colour code to easily recognise their diameter and a small apex diameter to create the correct implant preparation in both compact and soft bone. Milled and laser-marked working depth references. External irrigation drills. Diamond Like Carbon surface coating. The final drill, indicated in ONE surgical procedures, allows the preparation of the implant site to be completed with a countersink appropriate to the size of the implant.



### Optional drills

Optional drills make it possible to mark the spot where to start preparing the implant site (Marking drill) and to easily drill the cortical bone (Pointed starter drill). The Lindemann drill lets you correct the axis of the implant site after the use of the first surgical drill. The drill extension allows the length of contra-angle tools to be increased.



## Step drills stop set

Available in three different diameters to minimise the horizontal space occupied by the drill with the stop applied.



SET STOP DRILL CODE: OST10

### 28ST006

Drill stop  
 $\varnothing 2.0$  &  $\varnothing 2.4-2.8$  / L 6.0 mm

### 28ST008

Drill stop  
 $\varnothing 2.0$  &  $\varnothing 2.4-2.8$  / L 8.0 mm

### 28ST010

Drill stop  
 $\varnothing 2.0$  &  $\varnothing 2.4-2.8$  / L 10.0 mm

### 28ST011

Drill stop  
 $\varnothing 2.0$  &  $\varnothing 2.4-2.8$  / L 11.5 mm

### 28ST013

Drill stop  
 $\varnothing 2.0$  &  $\varnothing 2.4-2.8$  / L 13.0 mm

### 36ST006

Drill stop  
 $\varnothing 2.8-3.2$  &  $\varnothing 3.2-3.65$  / L 6.0 mm

### 36ST008

Drill stop  
 $\varnothing 2.8-3.2$  &  $\varnothing 3.2-3.65$  / L 8.0 mm

### 36ST010

Drill stop  
 $\varnothing 2.8-3.2$  &  $\varnothing 3.2-3.65$  / L 10.0 mm

### 36ST011

Drill stop  
 $\varnothing 2.8-3.2$  &  $\varnothing 3.2-3.65$  / L 11.5 mm

### 36ST013

Drill stop  
 $\varnothing 2.8-3.2$  &  $\varnothing 3.2-3.65$  / L 13.0 mm

### 45ST006

Drill stop  
 $\varnothing 3.65-4.1$  &  $\varnothing 4.1-4.5$  / L 6.0 mm

### 45ST008

Drill stop  
 $\varnothing 3.65-4.1$  &  $\varnothing 4.1-4.5$  / L 8.0 mm

### 45ST010

Drill stop  
 $\varnothing 3.65-4.1$  &  $\varnothing 4.1-4.5$  / L 10.0 mm


### 45ST011

Drill stop  
 $\varnothing 3.65-4.1$  &  $\varnothing 4.1-4.5$  / L 11.5 mm


### 45ST013

Drill stop  
 $\varnothing 3.65-4.1$  &  $\varnothing 4.1-4.5$  / L 13.0 mm


## DIAGRAM FOR POSITIONING AND REMOVING STOPS FROM DRILLS




1 **INSERT** the drill in the handpiece



2 **INSERT** the stop



3 **STOP** in position



4 **REMOVE** the drill stop

# COMPOSE YOUR KIT

## One Internal Implant Driver.


### ONE INTERNAL Implant Driver Set

The set includes the drivers suitable for implant connection with internal hexagon in mechanical and manual versions.

SET CODE: NDS10


#### SET COMPOSITION

<b>07NDD02</b>	Handpiece implant driver long L 11.0 mm
<b>02NDM10</b>	Implant manual driver short L 5.0 mm
<b>02NDM11</b>	Implant manual driver medium L 11.0 mm



### Implant drivers

The drivers for the ONE INTERNAL implant feature a retention system for engaging the implant, picking it up, transporting it and finalising the placement. Implant handpiece drivers feature the W&H hexagonal coupling system. They have laser reference markings indicating the transmucosal heights of the predefined prosthetic abutments (positioned at 2.0-3.5-4.5-5.5 mm). They have a reference that allows the hexagonal index of the implant to be phased appropriately (indicating the position of the vertex of the hexagon).



# COMPOSE YOUR KIT

## One Conical Implant Driver.


### ONE CONICAL Implant Driver Set

The set includes drivers suitable for the conical implant connection, available in mechanical and manual versions. See the ONE CONICAL catalogue section for all implant driver variants.

SET CODE: GDS10


#### SET COMPOSITION

<b>02GDM10</b>	Manual implant driver L 9.0 mm
<b>07GDD01</b>	Handpiece implant driver short L 9.0 mm
<b>07GDD02</b>	Handpiece implant driver long L 15.0 mm



### Implant drivers


The drivers feature a retention system to be able to engage the implant, pick it up, transport it and finalise the placement. All implant handpiece drivers feature the W&H hexagonal clamping system. They have laser reference markings indicating the transmucosal heights of the predefined prosthetic abutments (positioned at 2.0-3.5-4.5-5.5 mm). They also have a reference that allows the octagonal index of the implant to be appropriately phased (indicating the position of the vertex of the octagon).



# PROSTHETIC drivers.

## Prosthetic drivers

The ONE implant system has only one type of prosthesis driver (hexagon 1.2 mm) with which to screw and unscrew prosthetic components such as surgical cover screws, healing abutments, prosthetic component fixation screws and transfer abutments. Prosthesis handpiece drivers feature the W&H hexagonal coupling system.



**02EGM05**  
Manual prosthetic driver  
short  
L 3.0 mm




**02EGM10**  
Manual prosthetic driver  
medium  
L 10.0 mm




**02EGM15**  
Manual prosthetic driver  
long  
L 15.0 mm



**07-EG05**  
Handpiece prosthetic driver  
short  
L 5.0 mm




**07-EG10**  
Handpiece prosthetic driver  
medium  
L 10.0 mm




**07-EG15**  
Handpiece prosthetic driver  
long  
L 15.0 mm

## Dynamic driver

The ONE system features a prosthetic driver specifically for screwing at an angle to the implant axis. A useful solution for screw-retained prosthetic rehabilitations, this dynamic driver allows the driver to be angled up to 29° with respect to the implant axis.

 The dynamic driver is not compatible the standard fixing screws and must only be used with the relevant dynamic screws. Dynamic drivers are identified by a yellow colour-coded label.



**02ODM10**  
Dynamic driver short  
L 10.0 mm




**02ODM15**  
Dynamic driver long  
L 15.0 mm

## MUA and Overdenture driver

Specific drivers for straight MUA abutments and overdenture abutments. MUA drivers feature the W&H hexagonal clamping system.



**07DMU01**  
Handpiece MUA driver



**760CE**  
Driver OVD

# ACCESSORY TOOLS.

## Accessory tools

The ONE system features tools to facilitate specific surgical or prosthetic operations. The Direction Indicators are graduated and can therefore also be used to measure the depth of the osteotomy. The Depth Gauge is graduated and has a clasp tooth to check the position of the basal cortical in large sinus elevations. The torque ratchet has an adjustable torque control (10-70 Ncm) and can be used in both the surgical and prosthetic phases. The Ratchet Adapter allows manual use of all handpiece tools and can be used in conjunction with the torque ratchet. The Titanium Pliers allow contamination-free handling of the implant when needed.



**01-MP12**  
Double diameter  
depth probe  
(2 pieces)



**01-MP20**  
Double depth probe



**02-CT20**  
Torque ratchet  
10-70 Ncm



**02-AC50**  
Ratchet adapter



**01-PT10**  
Titanium forceps

## Bone tray


The Surgical Kit can be completed with a tray useful for handling autologous bone or bone substitutes.



**BOWL**  
Bone bowl

## X-ray slides


The X-ray slides allow for proper planning with CT, RX and OPT.



**01OLT10**  
ONE implant-sizing overlays







Implant line  
**packaging.**

---

# Implant line PACKAGING.

---

## Packaging

The packs are designed to be easy to use during surgery, providing clear identification of the Advan implant line, implant size, connection type or prosthesis type.

The dental implants are processed, cleaned and treated according to validated procedures. After decontamination, the dental implants are immediately stored. The implant screws will only be removed from storage in the Clean Room of the production plant. In this way, the surface oxide layer will form in a controlled atmosphere, not allowing contamination of the implant surface. After packing in the Clean Room, the implant screws are sent to the sterilisation process using  $\beta$ -rays.



## Visual identification of the implant

The pack cap allows the type of connection used to be recognised by the colour code:

- The system is packaged without a pre-assembled mounting device;
- The cover screw is included in the package;
- The self-centring geometry of the cover screw allows easy insertion and screwing;
- The cover screw is screwed in with Prosthesis Driver.


## Maximum Safety

- Decontamination by cold Argon plasma in Clean Room;
- Sterilisation by  $\beta$ -rays;
- Sterile barrier: PET-G blister sealed with TYVEK® foil;
- Internal packaging provided by the poly-carbonate vial;
- Cover screw positioned in the vial cap below a TYVEK® seal;
- Shaped blister to avoid shocks;
- Dental implant in contact only with titanium elements in order to avoid bimetallism phenomena or contamination from contact with plastic material.


Please refer to the ONE Surgical Manual for further information on the implant, packaging and surgical procedures.



## Label of non-sterile prosthetic packaging.



## Label of instruments packaging.



**Labels** of implant packs and sterile prosthetic components.

PRODUCT CODE — REF C3810 LOT XXFXXXX

COLOUR CODE — Ø=3,75 mm L=10,0 mm

BATCH NUMBER —

DIAMETER AND LENGTH OF THE IMPLANT —

PRODUCT FAMILY DESCRIPTION — ONE Conical

DATAMATRIX UDI CODE —

EXPIRY DATE (YEAR, MONTH, DAY) — YYYY-MM-DD

PRODUCT DESCRIPTION FOR DIFFERENT LANGUAGES —

PRODUCT CODE — REF 050MG12 LOT XXPXXXX

BATCH NUMBER —

ANTI-ROTATIONAL INDEX TYPE —

PRODUCT FAMILY DESCRIPTION — ONE Conical

DATAMATRIX UDI CODE —

EXPIRY DATE (YEAR, MONTH, DAY) — YYYY-MM-DD

PRODUCT DESCRIPTION FOR DIFFERENT LANGUAGES —

**Label** of implant blisters, sterile prosthetic components.

All the information on the implant blister label is replicated on the adhesive labels for the patient file.

PRODUCT CODE — REF C3810

BATCH NUMBER — LOT XXFXXXX

DIAMETER AND LENGTH OF THE IMPLANT — Ø=3,75 mm L=10,0 mm

EXPIRY DATE (YEAR, MONTH) — AAAAMM

COLOUR CODE —








**Label** of implant vial.

PRODUCT CODE — REF C3810

BATCH NUMBER — LOT XXFXXXX

## Indications for picking up the implant

### INSTRUCTIONS

1	<p>Choose implant type, length and diameter and take the blister out of the box.</p>	
2	<p>The container with the implant is sterile and contained in the blister. The label shows the product description and the batch number. Open the blister.</p>	
3	<p>Place the vial with the implant inside on a sterile surface.</p>	
4	<p>The cover screw is housed in the container cap under a heat-sealed Tyvek seal.</p>	
5	<p>Gently open the cap. (do not tear upwards)</p>	
6	<p>Connect the implant driver (choose the implant driver corresponding to the type of implant connection) and screw in at low speed (10-15 revolutions per minute).</p>	
7	<p>Remove the cover screw from the vial cap and screw it onto the implant using the Prosthetic Driver. (07-EG05; 07-EG10; 07-EG15; 02EGM10; 02EGM05; 02EGM15)</p>	



## Symbol description: packaging labels

### KEY



Manufacturer



Date of manufacture



Use-by date



Batch code



Catalogue number



Distributor



Sterilised using irradiation



Do not re-sterilise



Non-sterile



Do not use if the packaging is damaged and read the instructions for use



Single sterile barrier system with protective packaging inside



Single sterile barrier system with protective packaging outside



Keep away from sunlight



Keep dry



Do not re-use



Consult instructions for use or consult electronic instructions for use



Caution



Medical device



Unique device identifier



Not locking prosthetic component



Octagon locking prosthetic component



Hex locking prosthetic component



Multi packaging (the number reported in the symbol refers to the number of units in the packaging)



Advan products covered by the CE mark without the identification number fulfill the requirements of the Directive 93/42/CEE concerning medical devices and falls within Class I



Advan products covered by the CE mark fulfill the requirements of the Directive 93/42/CEE concerning medical devices and falls within Classes Im, Is, IIa, IIb

# IMPLANT CARD

## Advan.

### The patient Passport becomes SMART

The Implant Card is a real identity card that the customer will carry with him or her at all times, along with other personal documents.

This implant passport allows all practitioners to contact Advan directly for assistance and advice and to identify the components used, in order to intervene effectively.

This passport certifies the **originality** and **quality** of ADVAN implants and their components.




Please note that ADVAN Srl implant systems are manufactured in accordance with the ISO13485 quality management system certified by TÜV SÜD GmbH.

## Instructions for filling the passport out


### Implant Card

Front - Back for illustrative purposes only.



Information for filling out by the healthcare establishment/healthcare professional:

	Enter name or ID of the patient
	Name and address of the healthcare establishment/centre
	Date of Implantation
	Number of the tooth where the dental implant is positioned
	Name of the medical device
	Manufacturer
	Website with information for patients
	Batch number
	Product code
	Unique Device Identifier (UDI) in "machine-readable" format
<b>UDI-DI</b>	Unique Device Identifier (UDI) in "human-readable" format



**Magellan**

Guided surgery.

---

# MAGELLAN

## Computer Guided Surgery.

### Magellan

MAGELLAN is ADVAN's guided surgery system that allows all-round clinical case management, from software planning to final surgery.

The Magellan software application makes it possible to perform three-dimensional implant simulation directly on a personal computer. It also makes it possible to simulate the position of implants on two- and three-dimensional models, identify the mandibular canal, and draw overviews and sections of the bone model, while also offering the possibility of calculating bone density. By using Magellan, the dentist can plan implant-prosthetic surgery more safely, efficiently and quickly.

The surgical kit consists of a single box containing all the surgical and prosthetic tools necessary for performing the osteotomy, placing the implant and its prosthesis. Including surgical box made of Radel (non-toxic, autoclavable, metal-free).

The content of the Magellan guided surgery kit was developed by ADVAN to support the complete guided surgery workflow, from implant site preparation to final implant placement.

All tools and drills in the kit are used with a single Ø 4.6 mm guide sleeve without the use of stops or spanners to be assembled to the components at a later date, allowing surgery time to be reduced to a minimum.

The drills in the ONE Magellan kit allow the use of all implant diameters and lengths. In addition, in order to allow the correct insertion of ONE implants in a subcrestal position, special implant mounts, oversized in length by 1.5 mm, were made. The drills in the Magellan ONE kit allow the use of all implant diameters, with the exception of the 4.65 mm diameter, and all lengths.

## Magellan Guided Surgery Kit

The surgical kit consists of a single box containing all the surgical and prosthetic tools needed to perform the osteotomy and place the implant using a surgical template.

The kit contains the drivers necessary for **implant prosthesis**. Including surgical box made of Radel (non-toxic, autoclavable, metal-free).



OCG10	ONE GUIDED SURGERY KIT	QTY
02-AC50	Ratchet adapter	1
FS14	Lateral pin	2
FS15	Lateral pin drill	1
07-EG10	Prosthetic handpiece driver medium L 10.0 mm	1
07MCG10	Guided surgery mucotome	1
07DMG10	Guided surgery implant mount driver	1
02ECG10	Guided surgery mount extractor	1
28FC006	Surgical drill Ø 2.8 / L 6.0 mm	1
28FC008	Surgical drill Ø 2.8 / L 8.0 mm	1
28FC010	Surgical drill Ø 2.8 / L 10.0 mm	1
28FC011	Surgical drill Ø 2.8 / L 11.5 mm	1
28FC013	Surgical drill Ø 2.8 / L 13.0 mm	1
32FC006	Surgical drill Ø 3.2 / L 6.0 mm	1
32FC008	Surgical drill Ø 3.2 / L 8.0 mm	1
32FC010	Surgical drill Ø 3.2 / L 10.0 mm	1
32FC011	Surgical drill Ø 3.2 / L 11.5 mm	1
32FC013	Surgical drill Ø 3.2 / L 13.0 mm	1
36FC006	Surgical drill Ø 3.6 / L 6.0 mm	1
36FC008	Surgical drill Ø 3.6 / L 8.0 mm	1

OCG10	ONE GUIDED SURGERY KIT	QTY
36FC010	Surgical drill Ø 3.6 / L 10.0 mm	1
36FC011	Surgical drill Ø 3.6 / L 11.5 mm	1
36FC013	Surgical drill Ø 3.6 / L 13.0 mm	1
41FC006	Surgical drill Ø 4.1 / L 6.0 mm	1
41FC008	Surgical drill Ø 4.1 / L 8.0 mm	1
41FC010	Surgical drill Ø 4.1 / L 10.0 mm	1
41FC011	Surgical drill Ø 4.1 / L 11.5 mm	1
41FC013	Surgical drill Ø 4.1 / L 13.0 mm	1
07FPC10	Guided surgery crestal pin drill	1
01MCG01	One Conical - Guided surgery implant mount H 0.0 mm (2 pieces)	2
01MCN01	One Internal - Guided surgery implant mount H 0.0 mm (2 pieces)	2
07-FZ10	Guided surgery planisher	1
01PCG10	Guided surgery crestal pin	2
01VMG16	One conical - Guided surgery mount retaining screw (2 pieces)	1
01VMG18	One internal - Guided surgery mount retaining screw (2 pieces)	1
FN02	Guided surgery lateral pin sleeve (3 pieces)	1
01BCG45	Guided surgery sleeve (5 pieces)	1

## Magellan ONE drills

The drills for guided surgery in the ONE implant line have a single ring colour code to easily recognise their diameter. They have a small diameter apex to create the correct implant preparation with soft and spongy bone. The drills perform an osteotomy that can be considered congruous to the length of the chosen implant, the maximum increase in depth being only 0.4 mm (due to the conical apex of the drill as opposed to the convex apex of the implant). External irrigation drills.



28FC000	SET SURGICAL DRILL Ø 2.8 mm	QTY
<b>28FC006</b>	Surgical drill Ø 2.8 / L 6.0 mm	1
<b>28FC008</b>	Surgical drill Ø 2.8 / L 8.0 mm	1
<b>28FC010</b>	Surgical drill Ø 2.8 / L 10.0 mm	1
<b>28FC011</b>	Surgical drill Ø 2.8 / L 11.5 mm	1
<b>28FC013</b>	Surgical drill Ø 2.8 / L 13.0 mm	1



32FC000	SET SURGICAL DRILL Ø 3.2 mm	QTY
<b>32FC006</b>	Surgical drill Ø 3.2 / L 6.0 mm	1
<b>32FC008</b>	Surgical drill Ø 3.2 / L 8.0 mm	1
<b>32FC010</b>	Surgical drill Ø 3.2 / L 10.0 mm	1
<b>32FC011</b>	Surgical drill Ø 3.2 / L 11.5 mm	1
<b>32FC013</b>	Surgical drill Ø 3.2 / L 13.0 mm	1



36FC000	SET SURGICAL DRILL Ø 3.6 mm	QTY
<b>36FC006</b>	Surgical drill Ø 3.6 / L 6.0 mm	1
<b>36FC008</b>	Surgical drill Ø 3.6 / L 8.0 mm	1
<b>36FC010</b>	Surgical drill Ø 3.6 / L 10.0 mm	1
<b>36FC011</b>	Surgical drill Ø 3.6 / L 11.5 mm	1
<b>36FC013</b>	Surgical drill Ø 3.6 / L 13.0 mm	1



41FC000	SET SURGICAL DRILL Ø 4.1 mm	QTY
<b>41FC006</b>	Surgical drill Ø 4.1 / L 6.0 mm	1
<b>41FC008</b>	Surgical drill Ø 4.1 / L 8.0 mm	1
<b>41FC010</b>	Surgical drill Ø 4.1 / L 10.0 mm	1
<b>41FC011</b>	Surgical drill Ø 4.1 / L 11.5 mm	1
<b>41FC013</b>	Surgical drill Ø 4.1 / L 13.0 mm	1




## Magellan assembly tools and devices

Implant mount devices allow the correct positioning of the implant in the bone in both crestal and subcrestal positions, thanks to special implant mounts, oversized in length by 1.5 mm. Also included in the kit are the drivers for the mounting devices and the handy puller in case you have difficulty removing the mounting device from the system.


### ONE INTERNAL MAGELLAN

01KCG40	ONE INTERNAL TOOL SET	QTY
<b>07DMG10</b>	Guided surgery implant mount driver	1
<b>02ECG10</b>	Guided surgery mount extractor	1
<b>01MCN01</b>	One Internal - Guided surgery Implant mount H 0.0 mm (2 pieces)	2
<b>01VMG18</b>	One Internal - Guided surgery mount retaining screw (2 pieces)	2
<b>07EMG10</b>	Guided surgery manual mount driver	1




### ONE CONICAL MAGELLAN

01KCG30	ONE CONICAL TOOL SET	QTY
<b>01MCG01</b>	One Conical - Guided surgery Implant mount H 0.0 mm (2 pieces)	1
<b>01MCG02</b>	One Conical - Guided surgery Implant mount H 1.5 mm (2 pieces)	1
<b>01VMG16</b>	One Conical - Guided surgery mount retaining screw (2 pieces)	2
<b>07DMG10</b>	Guided surgery implant mount driver	1
<b>02ECG10</b>	Guided surgery mount extractor	1
<b>07EMG10</b>	Guided surgery manual mount driver	1



### ONE MAGELLAN

01KCG20	OPTIONAL DRILL AND PIN SET	QTY
<b>07MCG10</b>	Guided surgery mucotome	1
<b>07-FZ10</b>	Guided surgery planisher	1
<b>FS14</b>	Guided surgery lateral pin	2
<b>FS15</b>	Guided surgery lateral pin drill	1
<b>01PCG10</b>	Guided surgery crestal pin	2
<b>07FPC10</b>	Guided surgery crestal pin drill	1
<b>FN02</b>	Guided surgery lateral pin sleeve (3 pieces)	1



# LIFETIME WARRANTY

## Advan.

### Advan Warranty

Advan aims to provide the Professional with products that guarantee safety, effectiveness and enable professional efficiency. As part of the day-to-day development and improvement of its products, Advan reserves the right to make changes to products, tools, surgical techniques, prosthetic techniques and instructions for use if this constitutes an improvement for the product and an advantage for the Professional. Advan reserves the right to change prices, delivery times and terms and conditions of sale without prior notice. The products are not available in all countries.

Lifetime warranty: all Advan implants and prosthesis components are guaranteed against defects due to both the materials used and the manufacturing processes. This warranty only applies when Advan products are used in combination with original components and tools, used in accordance with the Instructions for Use (IFU). The warranty is only valid for the Dentist and the Dental Technician and cannot be passed on and challenged by the Patient. All methods of use of the products are indicated in the protocols of use and in the surgical manuals, which we recommend that you read carefully and that you scrupulously comply with the indications contained therein.

Advan replaces defective products depending on the product and the type of intervention required. The original buyer assumes all risks and liability arising from the use of the prod-

uct, whether used separately or together with products not distributed by Advan. The lifetime warranty is limited to the replacement of the defective product and will not reimburse the Dentist or Dental Technician or the Patient for any additional parts, tools, treatment costs or other costs and charges arising from or related to the replacement of the defective product. Treatment planning and the clinical application of the products are the sole responsibility of the dentist; Advan recommends attending intensive courses, obtaining a specialisation in implantology, and complying with the Instructions for Use (IFU) and the indications in the technical documentation available to users. Products may only be sold to Dentists and Dental Technicians in accordance with the requirements of the laws in force. The products in this catalogue may be covered by patents. Any replacement may only take place after filling out and sending to Advan the Return Form, available upon request from Customer Service.

The Return Form must be completed in its entirety, stamped and signed, accompanied by the possibly defective product and any other components installed by the dentist. In the event of implant failure, it is MANDATORY to enclose the X-ray documentation of the failed implant. Before sending any product, sterilise it.

**Important** - Read the instructions before use.

### Limits and Limitations

ADVAN DISCLAIMS ALL WARRANTIES, EXPRESS OR IMPLIED, AND EXCLUDES ALL LIABILITY FOR LOST PROFITS AND DIRECT OR INDIRECT DAMAGES, AS WELL AS COLLATERAL AND CONSEQUENTIAL DAMAGES, DIRECTLY OR INDIRECTLY RELATED TO ADVAN'S PRODUCTS, SERVICES OR INFORMATION.



LIVE YOUR LIFE, SMILE TO LIFE.

## BIBLIOGRAPHY

- 1. Baggi, L., Cappelloni, I., Di Girolamo, M., Maceri, F., & Vairo, G. (2008).** The influence of implant diameter and length on stress distribution of osseointegrated implants related to crestal bone geometry: A three-dimensional finite element analysis. *The Journal of Prosthetic Dentistry*, 100(6), 422–431.  
- [https://doi.org/10.1016/s0022-3913\(08\)60259-0](https://doi.org/10.1016/s0022-3913(08)60259-0)
- 2. Baggi, L., Di Girolamo, M., Mirisola, C., & Calcaterra, R. (2013b).** Microbiological evaluation of bacterial and mycotic seal in implant systems with different Implant-Abutment interfaces and closing torque values. *Implant Dentistry*, 22(4), 344–350.  
- <https://doi.org/10.1097/id.0b013e3182943062>
- 3. Biasotto, M., Scuor, N., Chiesa, R., Schimd, C., Di Lenarda, R., & Cigada, A. (2001).** Caratterizzazione di superfici di titanio per impianti. *Biomateriali*, 15(1), 5-13.
- 4. Calcaterra, R., M, D. G., Mirisola, C., & Baggi, L. (2016).** Effects of repeated screw tightening on implant abutment interfaces in terms of bacterial and yeast leakage in vitro: One-Time Abutment versus the Multiscrewing technique. *International Journal of Periodontics & Restorative Dentistry*, 36(2), 275–280.  
- <https://doi.org/10.11607/prd.2082>
- 5. Di Girolamo, M., Baggi, L., Pirelli, P., Pappalardo, L., Massei, G., Iaculli, F., Iezzi, G., Piattelli, A., & Calcaterra, R. (2017).** 3D X-ray microscopic analysis on a prosthetically loaded implant with platform-switching and conical connection: a case report. *Oral Implantology (Rome)*, 10(3), 241-246.  
- <https://doi.org/10.11138/orl/2017.10.3.241>
- 6. Di Girolamo, M., Calcaterra, R., Di Gianfilippo, R., Arcuri, C., & Baggi, L. (2017).** Bone level changes around platform switching and platform matching implants: a systematic review with meta-analysis. *PubMed*, 9(1), 1–10.  
- <https://doi.org/10.11138/orl/2016.9.1.001>
- 7. Hansson, S. (2000).** Implant-Abutment Interface: Biomechanical Study of Flat Top versus Conical. *Clinical Implant Dentistry and Related Research*, 2(1), 33–41.  
- <https://doi.org/10.1111/j.1708-8208.2000.tb00104.x>
- 8. Hansson, S. (2003).** A conical implant–abutment interface at the level of the marginal bone improves the distribution of stresses in the supporting bone. *Clinical Oral Implants Research*, 14(3), 286–293.  
- <https://doi.org/10.1034/j.1600-0501.2003.140306.x>
- 9. Meleo, D., Baggi, L., Di Girolamo, M., Di Carlo, F., Pecchi, R. & Bedini, R. (2012).** Fixture-abutment connection surface and micro-gap measurements by 3D micro-topographic technique analysis. *Ann Ist Super Sanità*, Vol. 48, No. 1: 53-58.  
- [https://doi.org/10.4415/ANN\\_12\\_01\\_09](https://doi.org/10.4415/ANN_12_01_09)
- 10. Rompen, E., Raepsaet, N., Domken, O., Touati, B., & Van Dooren, E. (2007).** Soft tissue stability at the facial aspect of gingivally converging abutments in the esthetic zone: A pilot clinical study. *The Journal of Prosthetic Dentistry*, 97(6), S119–S125.  
- [https://doi.org/10.1016/s0022-3913\(07\)60015-8](https://doi.org/10.1016/s0022-3913(07)60015-8)
- 11. Rocci, A., Calcaterra, R., Rocci, M., Rocci, C., Di Girolamo M., & Baggi, L. (2016).** Different performance of platform switching in equicrestal position implant: an histological study. *Oral Implantology (Rome)*, 9(1), 11-16.  
- <https://doi.org/10.11138/orl/2016.9.1.011>
- 12. Rocci, A., Calcaterra, R., Rocci, M., Rocci, C., & Baggi, L. (2017).** The influence of micro and macro-geometry in term of bone-implant interface in two implant systems: an histomorphometrical study. *PubMed*, 8(4), 87–95.  
- <https://doi.org/10.11138/orl/2015.8.4.087>
- 13. Small, P.N., & Tarnow, D.P. (2000).** Gingival recession around implants: a 1-year longitudinal prospective study. *The International Journal of Oral & Maxillofacial Implants*, 15(4), 527-32.  
- <https://doi.org/10.1097/00008505-200101000-00012>





SMILE TO LIFE  
WITH US.



**Advan S.r.l.**

Via Rosta della Maina, 2 - 33020 Amaro (UD)

+39 0433 096245

[www.advanimplantology.com](http://www.advanimplantology.com)